





<b>PHILIPS HEALTHCARE</b>	<b>Product Recycling Passport</b>	
-------------------------------	-----------------------------------	--

<b>Product name</b>	<b>Volcano CORE Series System</b>	
<b>Identification code(s)</b>	<b>400-0100.02</b>	
<b>Total weight (in Kg)</b>	<b>42.6</b>	
<b>Producer/ Manufacturer</b>	<b>Name company:</b>	<b>Philips Volcano</b>
	<b>Address:</b>	<b>2870 Kilgore Road Rancho Cordova, CA</b>
	<b>Zip code:</b>	<b>95670</b>
	<b>Country:</b>	<b>USA</b>

<b>Recycle Info</b>	<b>Items:</b>	<b>Location</b>
<b>Special attention</b>  	<b>N/A</b>	

  <b>Fluids / Gases</b>	<b>Items:</b>	<b>Location</b>
	<b>N/A</b>	

<b>Batteries To be Removed</b>  	<b>Type:</b>	<b>Location</b>
	<b>CR2032 Button Cell Lithium Battery 3.0V</b>	<b>Motherboard</b>

<p><b>Hazardous</b></p>  <p><b>To be Removed</b></p>	<b>Substances:</b>	<b>Location</b>
	may contain lead	Circuit board
	may contain hexavalent chromium	Screws, fasteners





Figure 5: Typical Volcano system



Coin type motherboard battery to be removed



### Recycling Passport Revision History

Rev	Description	Date
A	CORE System WEEE Passport	4/7/16