





|                               |                                   |  |
|-------------------------------|-----------------------------------|--|
| <b>PHILIPS<br/>HEALTHCARE</b> | <b>Product Recycling Passport</b> |  |
|-------------------------------|-----------------------------------|--|



|                                   |                          |   |
|-----------------------------------|--------------------------|---|
| <b>Product name</b>               | <b>Volcano HCB</b>       |   |
| <b>Identification code(s)</b>     | <b>409-0100.01</b>       |   |
| <b>Total weight (in Kg)</b>       | <b>2</b>                 |   |
| <b>Producer/<br/>Manufacturer</b> | <b>Name<br/>company:</b> | <b>Philips Volcano</b>                      |
|                                   | <b>Address:</b>          | <b>2870 Kilgore Road Rancho Cordova, CA</b> |
|                                   | <b>Zip code:</b>         | <b>95670</b>                                |
|                                   | <b>Country:</b>          | <b>USA</b>                                  |

| <b>Recycle Info</b>  | <b>Items:</b> | <b>Location</b> |
|--|---------------|-----------------|
| <b>Special attention</b><br><br> | <b>N/A</b>    |                 |
|  |               |                 |
|  |               |                 |
|  |               |                 |
|  |               |                 |

|  |               |                 |
|--|---------------|-----------------|
| <br><br><b>Fluids / Gases</b> | <b>Items:</b> | <b>Location</b> |
|  | <b>N/A</b>    |                 |
|  |               |                 |
|  |               |                 |
|  |               |                 |

|   |              |                 |
|---|--------------|-----------------|
| <br><br><b>Batteries</b><br><br><b>To be Removed</b> | <b>Type:</b> | <b>Location</b> |
|   | <b>N/A</b>   |                 |
|   |              |                 |
|   |              |                 |
|   |              |                 |

|   |                                 |                   |
|---|---------------------------------|-------------------|
| <p>Hazardous</p>  <p>To be Removed</p> | Substances:                     | Location          |
|   | may contain lead                | Circuit board     |
|   | may contain hexavalent chromium | Screws, fasteners |
|   |                                 |                   |
|   |                                 |                   |
|   |                                 |                   |



### Recycling Passport Revision History

| Rev | Description          | Date   |
|-----|----------------------|--------|
| A   | FM-Hub WEEE Passport | 4/7/16 |