





PHILIPS HEALTHCARE	Product Recycling Passport	
-------------------------------	-----------------------------------	--

Product name	Volcano s5i Imaging System Workstation	
Identification code(s)	806300-001	
Total weight (in Kg)	70	
Producer/ Manufacturer	Name company:	Philips Volcano
	Address:	2870 Kilgore Road Rancho Cordova, CA
	Zip code:	95670
	Country:	USA

Recycle Info	Items:	Location
Special attention 	N/A	

 Fluids / Gases	Items:	Location
	N/A	

 Batteries To be Removed	Type:	Location
	CR2032 Button Cell Lithium Battery 3.0V	Motherboard

Hazardous  To be Removed	Substances:	Location
		may contain lead
	may contain hexavalent chromium	Screws, fasteners




Figure 5: Typical Volcano system



**Coin type
motherboard battery
to be removed**



Recycling Passport Revision History

Rev	Description	Date
A	S5i WEEE Passport	4/7/16